

The Sercon μ Carbs

The next generation acid delivery equilibration system for carbonate and DIC analysis

Key features include:

- Dual needle system with separate needles for acid delivery and headspace flush / sample transfer which optimise flow rates and minimise blockages
- Temperature programmable heated sample tray optimises reaction rates ($\pm 0.1^\circ\text{C}$)
- Sample tray holds up to 170 6 ml or 12 ml vials which can be analysed in a single, unattended run. Alternative trays can be provided for larger vials
- Phosphoric acid reservoir held within the heated sample tray ensures smooth acid delivery and minimises acid crystallisation
- Acid resistant lids manufactured from Ertalyle on the sample trays minimises accidental acid corrosion
- The hyphenated μ Carbs-HS2022 system generates high precision isotopic data from carbonate and dissolved inorganic carbon (DIC) samples. The high sensitivity source in the HS2022 ensures excellent precision at extremely small sample sizes
- The flexible system also allows the stable isotopes of water to be analysed via water equilibration and gas analysis such as CO_2 in breath or atmospheric samples

Sample measurement:


- $\delta^{18}\text{O}$ and $\delta^{13}\text{C}$ in CO_2 in carbonates
- $\delta^{13}\text{C}$ in dissolved inorganic carbon (DIC)
- $\delta^2\text{H}$ and $\delta^{18}\text{O}$ in water and other aqueous samples
- $\delta^{13}\text{C}$ on breath and other gaseous samples, including atmospheric samples



The Sercon μ Carbs. The next generation acid delivery equilibration system for carbonate and DIC analysis

Carbonate and DIC analysis


The sample vial is placed inside the heated sample tray set to the required temperature. The headspace of the vial is flushed with helium, then a second needle delivers phosphoric acid, in order that CO_2 be generated via the reaction below. After the equilibration time has elapsed, the headspace is transferred to the IRMS for isotopic analysis.




Water and other equilibration measurements

The sample vial is placed inside the heated sample tray set to the required temperature. The headspace of the vial is flushed with a mixture of helium and CO_2 or H_2 (depending upon which isotope is to be measured), in order that the reaction below can take place inside the vial. After the equilibration time has elapsed, the headspace is transferred to the IRMS for isotopic analysis

Oxygen isotope equilibration



Hydrogen isotope equilibration



Gas analysis

No sample pre-treatment is required and the sample can be measured directly. The sample is simply transferred from the vial to the IRMS via the sampling needle.



Autosampler

Sercon produce the most elegant autosamplers using advanced engineering technologies. The design is robust, reliable, easy to operate and maintain.

Service and Support

At Sercon we pride ourselves on the support available to customers, and consistently receive good feedback on the service we provide. We offer support via telephone, email and remote log on. If an engineer visit is necessary we provide rapid on-site response from our team of specialist, experienced engineers.

All users receive training as part of the installation programme. We can also provide further training on specific applications and tailor your course to your analytical needs.

Power and Gas Requirements

Power	100-240 VAC
Helium	99.998%

Specifications

Sample type	Sample size	Isotope	Precision (1 σ)
Carbonate	500 μg	$\delta^{13}\text{C}$	0.08 ‰
		$\delta^{18}\text{O}$	0.08 ‰
	100 μg	$\delta^{13}\text{C}$	0.08 ‰
		$\delta^{18}\text{O}$	0.1 ‰
	50 μg	$\delta^{13}\text{C}$	0.1 ‰
		$\delta^{18}\text{O}$	0.12 ‰
DIC	15 ppm	$\delta^{13}\text{C}$	0.15 ‰
Water	200 μl	$\delta^{18}\text{O}$	0.1 ‰
		$\delta^2\text{H}$	2 ‰
Atmospheric CO_2	12 ml	$\delta^{13}\text{C}$	0.1 ‰
		$\delta^{18}\text{O}$	0.2 ‰
Breath CO_2	12 ml	$\delta^{13}\text{C}$	0.1 ‰

Tel: +44 (0)1270 580008

Email: sales@sercongroup.com

www.sercongroup.com

June 2019

