



sercon
innovators in isotopes

ABCA2 Upgrade



Sercon System Controller Electronics Upgrade

New electronics – high performance, state of the art electronics, in a standalone unit which can be placed alongside the instrument

New PC and software included - enabling compatibility with all major data suites and programs, including all versions of MS Windows

Optional source and collector re-builds - for ultimate system performance.

XYZ Autosampler Upgrade

Robust and reliable XYZ sampler with five 48 place racks for 12ml vials meaning a total of 240 vials can be loaded at a time.

Optional Bar Code reader for secure patient identification and GLP compliance. With the bar code reader installed samples can be reloaded during a run, enabling a high number of consecutive analyses.



Pumping System Upgrade

Replacement of old Pfeiffer or Edwards turbopump with the latest Edwards nEXT85. Utilising existing gauge* and rotary pump, therefore there is no change to the instrument footprint.

*Gauge can be upgraded dependent upon lifetime, please discuss with the Sercon team

